



# chromoscreen®

## ChromoScreen® pack Content and Storage

---

In the ChromoScreen® pack you will find:



5 ChromoScreen slides




Compact disc with:  
ChromoScreen users' manual;  
file.gal; file.qc



50 Coverslips



ChromoScreen Labelling reagents



ChromoScreen Layout



ChromoScreen References

Store the Slides at Room Temperature in a dry environment. Store the ChromoScreen® Labelling reagents, COT Human DNA and References at -30 to -15°C in a nofrost- freezer. Use ChromoScreen® microarrays and all reagents provided before their expiry date.

More info and contacts at  
[www.chromoscreen.com](http://www.chromoscreen.com)

For research use only.  
Not for use in diagnostic procedures.

**IBSA Institut Biochimique SA**

Via del Piano 29, CH-6915 Pambio-Noranco, Lugano, Switzerland  
[www.ibsa-international.com](http://www.ibsa-international.com)

